
LIPEX[®] Flow

LIPEX[®] Flow Thermobarrel Extruder packages

OVERVIEW

LIPEX[®] Flow Thermobarrel Extruders are the industry-leading bench-top extruders for R&D and cGMP manufacturing of liposomal formulations in both academia and industry.

The innovative and patent-pending design maximizes the effective filtration area, resulting in lower extrusion pressures and high-flow processes. This translates to increased throughput, a wider formulation application range, faster processing times, and minimized process risks



*EFA = Effective Filtration Area **MAOP = Maximum Allowable Operating Pressure

The units are designed to produce a homogenous population of large unilamellar vesicles from a non-homogenous population of multi-lamellar vesicles. The unilamellar vesicles are formed by utilizing a constant pressure force of between 100 and 2400 psi to force the vesicles through filters of a predefined pore size. These units have a maximum operating pressure and temperature of 2450 psi and 80 °C, respectively. The side walls of the thermobarrel and filter support are jacketed to provide efficient temperature control of the extruder for optimal size reduction results.

PACKAGE CONTENTS

- (1x) LIPEX[®] Flow Thermobarrel Extruder
- (1x) Package of polyester drain discs (100/pack)
- (1x) Package of polycarbonate filters with 100 nm pore size (100/pack)
- (1x) Spare O-ring set
- (1x) High pressure nitrogen line
- (1x) Turnover package USB
 - Operating manual
 - Declaration of conformity
 - Certificate of inspection & testing



TECHNICAL SPECIFICATIONS

LIPEX [®] Flow Model	10 mL	100 mL	1000mL
MAOP*			
• psig	2450	2450	2450
• bar	170	170	170
MAOT (°C)**	80	80	80
Min. Extrusion volume	1 mL	20 mL	100 mL
Max. Extrusion volume	10 mL	100 mL	1000 mL
Weight (kg)	3	7	25
Material		•	•

• Body

• O-Ring

Surface finish

316 L Stainless steel EPDM

Mechanically finished



*MAOP = Maximum Allowable Operating Pressure **MAOT = Maximum Allowable Operating Temperature

REQUIRED ANCILLARY EQUIPMENT/SUPPLIES

The following is a list of commonly used ancillary equipment and materials required for efficient operation of the extruder.

- polyester drain disc (included)
- heated re-circulator
- high pressure gas cylinder (nitrogen)
- polycarbonate filter (included) gas regulator with a minimum delivery pressure of 2500 psi

This information and all further technical advice are based on our present knowledge and experience. However, it implies no liability or other legal responsibility on our part, including with regard to existing third party intellectual property rights, especially patent rights. In particular, no warranty, whether express or implied, or guarantee of product properties in the legal sense is intended or implied. We reserve the right to make any changes according to technological progress or further developments. The customer is not released from the obligation to conduct careful inspection and testing of incoming goods. Performance of the product described herein should be verified by testing, which should be carried out only by qualified experts in the sole responsibility of a customer. Reference to trade names used by other companies is neither a recommendation, nor does it imply that similar products could not be used.

LIPEX® Flow - reg. trademark of Evonik Industries AG and its subsidiaries

Evonik Operations GmbH Health Care Business Line Drug Delivery & Medical Device Solutions

healthcare@evonik.com www.evonik.com/healthcare

